



# Adding & Managing Users



# **A guide to adding & managing users in Real-Time Online**

This guide will take you through the steps of adding new users in Real-Time Online and managing their details, security roles and permissions within the system.


# Step 1

Once you've logged in to Real-Time Online, **click on 'Admin'** in the side menu and select **'Manage Users'** from the drop-down menu.

# Step 2


You'll be taken to a page that displays your current users.

If you're adding a new user, **click one of the 'Add User' buttons**. For ease, you'll find one of these at the top of the page and one of them below your table of users.



The screenshot shows the Invisible Systems dashboard. The left sidebar menu is open, and the 'Admin' option is selected, with a red arrow pointing to it. The 'Manage Users' option is highlighted with a red circle. The main content area displays 'Realtime Online - Medical organisations' with a table of data.

Organisations	Total Sensors	...With No Notifications
Clinics	7	2
Hospitals	26	6



The screenshot shows the 'User Management' page. The 'Add User' button is highlighted with a red circle. The page displays a table of users with columns for Enabled, User Name, Email Address, Mobile No., Phone No., Plain Text Email, and Text to Speech No.

Enabled	User Name*	Email Address*	Mobile No.	Phone No.	Plain Text Email	Text to Speech No.
Yes	Medical Demo	*****		*****	*****	
Yes	Healthcare Show 2022	*****				
Yes	Medical Demo Admin	support@invisible-systems.com				
Yes	George Anderson	*****				

# Step 3

From here, you'll be able to enter the new **user's name and email address**, plus any other relevant information.

## You can:

- Enter their **mobile phone number** if you want the user to receive SMS alerts when sensor alarms are triggered
- Add the user to their **department/organisation**, so they can only view specific sensors and data
- Click 'Application' to **choose the default dashboard** that the user will see when they log in to Real-Time Online
- Check the box 'Organisation Admin' if you want the user to have **access to any new systems** that are added to their organisation/department

Clinic 1

Add

Active

**User Name**

**Email Address**

We'll never share your email with anyone else.

**Mobile**

**Phone Number**

**Text to Speech**

**Plain text email address (Pagers & DECT phones)**

**Application**

**Organisation**

Always notify for alarm suppression periods

Can Change Password

**Password Change Interval**

**Password Expiry Date**


Cancel Save changes

Check the box 'Can change password' if you want the user to be able to **generate or change their own password**.

For most new users, this box will need to be ticked to allow them to generate a password for the first time.

You can also **set password intervals and expiry dates**. These are useful if you want users to change their passwords after a certain period of time for cyber security reasons, or if you're setting up a temporary user.

Click 'Save changes'.



The screenshot shows a user profile form titled "Add" for "Clinic 1". The form includes several sections: "Active" (checked), "User Name" (text input), "Email Address" (text input with a note "We'll never share your email with anyone else."), "Mobile" (text input with a dropdown for country code), "Phone Number" (text input with a dropdown for country code), "Text to Speech" (text input with a dropdown for country code), "Plain text email address (Pagers & DECT phones)" (text input), "Application" (dropdown menu set to "Mixed Dashboard"), and "Organisation" (dropdown menu set to "Clinics"). At the bottom, there is a checkbox for "Always notify for alarm suppression periods" (unchecked) and a section highlighted with a pink border containing "Can Change Password" (checked), "Password Change Interval" (dropdown set to "Never"), and "Password Expiry Date" (text input set to "Password Expires"). At the bottom right of the form are "Cancel" and "Save changes" buttons.

# Step 4

You will then see this permissions screen, where you can select the roles that the new user will have.

Leave the boxes unticked if you want your user to be **'read-only'**, i.e. they will only be able to view data - not edit anything.

A **'System User'** can access data, add new audits and acknowledge notifications.

A **'System Admin'** can make changes to the system, such as adding/editing alarms and sensors, changing the names of sensors, suppressing alarms, etc.

An **'Organisation Admin'** can manage organisation-wide activities, like adding new systems and organisations.

Click **'Save changes'**.

	System User	System Admin	Organisation Admin
Select all	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Core Features	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Conditions Monitoring	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
View Legacy Pages	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

	System User	System Admin	Organisation Admin
Select all	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Core Features	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Conditions Monitoring	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
View Legacy Pages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>


	System User	System Admin	Organisation Admin
Select all	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Core Features	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Conditions Monitoring	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
View Legacy Pages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

# Step 5

A User Site Access window will appear.

**Choose which systems the user is allowed to access** by checking the box next to each one. You can use the search bar to search for specific systems.

**Click 'Save changes'.**




# Step 6

You will see this window asking if you want to send a **password reset email** to the new user. **Click 'Ok'** to send this now.

The email will be sent to the user's inbox, where they can then choose a password and log in to Real-Time Online.

If you don't want to send a password reset email to the new user, simply **click 'Cancel'**.



# Step 7



Once you have added your new user, you will be taken back to the screen displaying all of your current users. Your new user will now have been added to this table.

If you need to edit a user's details, permissions or site access at any point, **click on the cog button** next to that user under 'Actions'.

To find a specific user in a system, type their name into the search bar.

# Step 8

If you ever want to **deactivate a user**, you can do so by clicking 'Edit User', then **deselecting the 'Active' checkbox** at the top of the window. You can reactivate them at any time by reticking the box.





If you require any additional training, please contact  
[help@invisible-systems.com](mailto:help@invisible-systems.com)

